

Production of lentiviruses in HEK293T cells 2nd Generation lentiviral system

1. Background and purpose of the procedure

Lentiviral transduction is an effective method for creating a stable cell line with a DNA cassette of interest integrated into its genomic DNA. The second-generation system has one packaging plasmid which includes all the important packaging components: Gag, Pol, Rev, and Tat. To produce virus, it requires a single packaging plasmid, an envelope plasmid, and a transfer vector (i.e., with the gene of interest to be expressed).

2. Materials and equipment

- HEK293T cells (ATCC, CRL-3216)
- Complete cell culture medium (DMEM, Gibco, 11995073) with FBS (FB Essence, VWR, 10803-034) and antibiotics (Penicillin/Streptomycin, Corning, 30-002-CI)
- D-PBS, without Ca⁺⁺ and Mg⁺⁺ (Gibco, 14190144)
- Opti-MEM Serum-free media (Gibco, 11995)
- Lentiviral packaging vector psPAX2 (Addgene, 12260)
- VSV-G envelope vector pMD2.G (Addgene, 12259)
- pLV transfer vector containing the gene of interest
- T75 TC-treated Flask
- 15 mL and 50 mL polypropylene tubes
- Cell counter/Hemocytometer

3. Method

1. Day 1

- a. For each plasmid to be transfected, plate 3×10^6 HEK293T cells in 10-15 mL of media in a T75 TC-treated flask
- b. Incubate the cells at 37°C, 5% CO₂ overnight

2. Day 2

- a. Perform the transfection in the late afternoon, transfection mix should only be incubated with the cells for 15-18 h maximum
- b. In a 1.5 mL microcentrifuge tube, mix:
 - i. 1 mL of Opti-MEM
 - ii. 10 µg of transfer vector (pLV)
 - iii. 10 µg of packaging vector (psPAX2)
 - iv. 5 µg of envelope vector (pMD2.G)
 - v. 75 µL of Lipofectamine reagent (3 µL:1 µg of DNA) [*yellow cap tube*]
 - vi. 50 µL of Lipofectamine 3000 (2 µL:1 µg of DNA) [*red cap tube*]
- c. Vortex for 10 sec
- d. Incubate the mix for 15-30 min at room temperature

- e. Gently add the DNA:Lipofectamine 3000 mix to the cells being careful to not detach the cells from the flask
- f. Swirl the flask gently to disperse the mixture evenly
- g. Incubate the cells for 15-18 h at 37°C, 5% CO₂

3. Day 3

- a. In the morning, change the media to remove the transfection reagent. Replace with 20 mL of fresh complete DMEM. Be careful to not detach the cells from the flask
- b. Incubate the cells for 48 h at 37°C, 5% CO₂

4. Day 5

- a. Harvest media containing the lentiviral particles and transfer to a 50 mL polypropylene tube
- b. Spin media at 1,500 rpm for 5 min any cells or cells debris that were collected during harvesting
- c. Viruses may be stored at 4°C for short-term storage (1 day to 1 week) or frozen at -80°C for long term-storage (up to 1 year, aliquots into smaller tubes to minimize freeze/thaw cycles)
- d. Add 10 mL of fresh complete DMEM. Be careful to not detach the cells from the flask [optional]
- e. Incubate the cells for 48 h at 37°C, 5% CO₂ [optional]

5. Day 6 [optional]

- a. Harvest media containing the lentiviral particles and transfer to a 15 mL polypropylene tube
- b. Spin media at 1,500 rpm for 5 min any cells or cells debris that were collected during harvesting
- c. Viruses may be combined with the viruses harvested on Day 5 and stored at 4°C for short-term storage (1 day to 1 week) or frozen at -80°C for long term-storage (up to 1 year, aliquots into smaller tubes to minimize freeze/thaw cycles)

4. References

Chu C, Xin A, Zhou Y, Zhang Y. A simple protocol for producing high-titer lentivirus. Acta Biochim Biophys Sin (Shanghai). 2013 Dec;45(12):1079-82. doi: [10.1093/abbs/gmt112](https://doi.org/10.1093/abbs/gmt112). Epub 2013 Oct 10. PMID: [24113089](https://pubmed.ncbi.nlm.nih.gov/24113089/).

[Producing lentivirus in HEK293T cells using a 2nd Generation lentiviral system](#), from David P. Turner, Medical University of South Carolina, Hollings Cancer Center

5. Revision history

Revision #	Date	Prepared by
1.0	2021-10-07	Elie Besserer-Offroy
Summary of modifications		
Initial version of the protocol		